

## Technical data sheet

Ref: FT LM-H1147

Page : 1/2

Version date: 24/07/2024

## **Minimum Essential Media**

w/ Hanks' Salts w/o L-Glutamine w/ 25mM Hepes

**CAT N°: LM-H1147** 

**Storage conditions**: +2°C to +8°C, protected from light

Shelf life: 24 months

**Composition**: Displayed on website; also available on request

Colour: Red, orange solution

**pH**: 7.2 ± 0.3

**Osmolality**: 300 mOsm/kg  $\pm$  10 %

**Endotoxin**: < 1 EU/ml

## **Sterility tests**:

- Bacteria in aerobic and anaerobic conditions

- Fungi and yeasts

**Cell Growth test:** Medium tested for the ability to support L929 or MRC5 cell growth

Other tests: Not applicable

#### Recommended use:

- Respect storage conditions of the product
- Do not use the product after its expiry date
- Store product in an area protected from light (not necessary for saline solutions).
- Manipulate the product in aseptic conditions (e.g.: under laminar air flow)
- Wear clothes adapted to the manipulation of the product to avoid contamination (e.g.: gloves, mask, hygiene cap, overall...)

The product is intended to be used in vitro for research or further manufacturing only and not for use as an Active Pharmaceutical Ingredient or food or animal feed.

### **Application:**

Minimum Essential Medium (MEM) is a modification of Eagle's earlier medium Basal Medium Eagle (BME), containing higher concentrations of the essential nutrients. Earle's salts have been replaced with Hanks' salts. The bicarbonate concentration is substantially lower, but this medium contains Hepes to have a good buffering capacity.



# Technical data sheet

Ref: FT LM-H1147

Page : 2/2

Version date : 24/07/2024

### **Uses**:

Supplements, such as antibiotics, should be added as sterile supplements to the medium. Storage conditions and shelf-life of supplemented product will be affected by the nature of the supplements. Sterile serum should not be refiltered before or after being added to sterile medium because growth promoting capacity may be reduced upon re-filtration. Add 10 ml/l of L-Glutamine 200 mM (BioSera CAT N $^{\circ}$ : XC-T1715) before using this medium. BioSera also recommend the addition of 10% of Fetal Bovine Serum in the medium

### **Signs of Deterioration:**

Medium should be clear and free of particulate and flocculent material. Do not use if medium is cloudy or contains precipitate.

Other evidence of deterioration may include colour change or degradation of physical or performance characteristics.

**Remarks:** Not applicable