

Iron supplementation treatment on sera

CAT N°: XX-xxxxR

For iron supplemented product the **letter R** is added at the end of the regular CAT n° serum whatever the specie.

E.g.: FB-1001 => FB-1001R / CA-115 => CA-115R

Analyses:

Tests and results are recorded on CoA.

Treatments:

A 1% Ferric citrate sterile solution is added to the sterile serum at 0.246% (v/v).

Then the serum is well mixed.

Transferrin is normally 20 to 40% bound to iron.

Iron concentration in treated product must be > 500 µg/100ml.

Remarks:

All other parameters than the ones mentioned above remain equal and can be find in each regular (non-treated) serum TDS.