

COMPOSITION

Trypsin		CAS number	LM-T1703 Liquid mg/l	LM-T1704 Liquid mg/l	LM-T1705 Liquid mg/l	LM-T1720 Liquid mg/l	LM-T1721 Liquid mg/l
Inorganic Salts	EDTA Disodium Salt Dihydrate	6381-92-6	/	/	221,4	380	200
	Potassium Chloride	7447-40-7	200	200	400	400	400
	Potassium Phosphate Monobasic Anhydrous	7778-77-0	200	200	/	60	60
	Sodium Bicarbonate	144-55-8	/	/	580	350	350
	Sodium Chloride	7647-14-5	8000	8000	8000	8000	8000
	Sodium Phosphate Dibasic Anhydrous	7558-79-4	1150	1150	/	48	48
O.C.*	D-Glucose Anhydrous	50-99-7	1000	1000	1000	1000	1000
	Phenol Red Sodium Salt	34487-61-1	2	/	2	10	10
	Trypsin 1:250	9002-07-7	2500	2500	500	2500	2500

Trypsin		CAS number	XC-T1716 Liquid mg/l	XC-T1726 Liquid mg/l	LM-T1706 Liquid mg/l	XC-T1717 Liquid mg/l	PM-T1717 Powder mg/l	LM-T1701 Liquid mg/l
Inorganic Salts	EDTA Disodium Salt Dihydrate	6381-92-6	/	/	254,8	2214	2548	372
	Potassium Chloride	7447-40-7	200	400	200	4000	/	200
	Potassium Phosphate Monobasic Anhydrous	7778-77-0	200	60	200	/	/	200
	Sodium Bicarbonate	144-55-8	/	350	/	/	/	/
	Sodium Chloride	7647-14-5	8000	8000	7950	8000	8040	7950
	Sodium Phosphate Dibasic Anhydrous	7558-79-4	1150	47,88	1150	/	/	1150
O.C.*	D-Glucose Anhydrous	50-99-7	1000	1000	1000	10000	/	1000
	Phenol Red Sodium Salt	34487-61-1	/	/	/	/	/	/
	Trypsin 1:250	9002-07-7	25000	25000	500	5000	500	/
	Recombinant Trypsin 2500 USP	9002-07-7	/	/	/	/	/	50

* Other Components